



SUCCEEDER®

Coagulation Reagents



ISO

ISO13485:2016



CE

No. Q5 076148 0008 Rev.03



Coagulation Reagents

1 Good Linearity And Sensitivity

- Good sensitivity of PT for factor deficiency, liver diseases and anticoagulation therapy monitoring.
- APTT reagent (lyophilized) main component is silica, with good linearity and high sensitivity.
- Wide detection linearity of FIB up to 0.8 – 8 g/L.


2 Good Stability After Reconstituted

- 2 – 8°C, stable for 10 days after reconstituted.

3 Reliable Quality

- GMP Factory with 100.000 class.
- ISO 13485 certificate.

4 Correlation



5 Specification

Product	Packing (per unit)	Unit
Prothrombin Time Kit (PT)-Liquid	12x2 ml/vial	Kit
Prothrombin Time Kit (PT)-Lyophilized	6x2 ml/vial	Kit
Prothrombin Time Kit (PT)-Liquid	12x4 ml/vial	Kit
Prothrombin Time Kit (PT)-Lyophilized	6x4 ml/vial	Kit
Activated Partial Thromboplastin Time Kit (APTT)-Liquid	6x2 ml/vial	Kit
Activated Partial Thromboplastin Time Kit (APTT)-Lyophilized	6x2 ml/vial	Kit
Activated Partial Thromboplastin Time Kit (APTT)-Liquid	6x4 ml/vial	Kit
Activated Partial Thromboplastin Time Kit (APTT)-Lyophilized	6x4 ml/vial	Kit
Fibrinogen Reagent Kit (FIB)-Liquid	6x2 ml/vial	Kit
Fibrinogen Reagent Kit (FIB)-Lyophilized	6x2 ml/vial	Kit
Fibrinogen Reagent Kit (FIB)-Liquid	6x4 ml/vial	Kit
Fibrinogen Reagent Kit (FIB)-Lyophilized	6x4 ml/vial	Kit
Thrombin Time Kit (TT)-Liquid	12x2 ml/vial	Kit
Thrombin Time Kit (TT)-Lyophilized	12x2 ml/vial	Kit
Thrombin Time Kit (TT)-Liquid	12x4 ml/vial	Kit
Thrombin Time Kit (TT)-Lyophilized	12x4 ml/vial	Kit
Coagulation Control Level I Kit	12x1ml/vial	Kit
Coagulation Control Level II Kit	12x1ml/vial	Kit
D-Dimer Kit (DD)-Liquid	R1:4x5.5 ml/vial R2:4x2 ml/vial	Kit
D-Dimer Kit (DD)-Liquid	R1:6x11 ml/vial R2:6x4 ml/vial	Kit
Fibrinogen Degradation Product Kit (FDP) -Liquid	R1:4x6 ml/vial R2:4x2 ml/vial	Kit
D-Dimer Control Level I Kit	6x0.5 ml/vial	Kit
D-Dimer Control Level II Kit	6x0.5 ml/vial	Kit

Thrombosis
Hemostasis



BEIJING SUCCEEDER TECHNOLOGY INC.

8/F, Tower 1, No.7 Kexueyuan Road, Changping District, Beijing 102206
 Tel: +86-10-53855569 Fax: +86-10-53855570
 E-mail:sales@succeeder.com.cn Website: www.succeeder.com